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# Manual

## Power Supply for mercury lamps (ST-Series)

DURATEC Analysentechnik GmbH  
Rheinauer Strasse 4  
DE-68766 Hockenheim  
GERMANY

Tel. +49-6205 / 9450-0  
Fax. +49-6205 / 9450-33  
eMail [info@duratec.com](mailto:info@duratec.com)  
Internet [www.duratec.com](http://www.duratec.com)

## 1 Packing List

1 pc	Power Supply
1 pc	Manual
1 pc	Power Cord
1pc	Lamp Cord

## 2 Specifications

Lamp current	: approx. 0 ... 1.0 Ampere
Ignition Voltage	: > 600 Volt
Supply	: 100/120/220/240VAC; 50/60 Hz
Fuses	: Lamp current 1.6 A slow blow (internal)
	AC input 1.60A @220/240V
	3.15A @100/120V
Dimensions	: W 22 cm x D 15 cm x H 35 cm
Weight	: 8 kg


## 3 Description

The power supply is designed for the use of mercury lamps with an anode current of 1 amp and a tube drop voltage of 120 volts or less. The specifications of this device are adjusted to this lamp.

The power supply delivers the ignition voltage and lamp current. Therefore the operation confines itself to starting the device and adjusting the lamp-current.

## 4 Parts identification and function


### 4.1 Front Panel

**Fig. 1 Front Panel**

Current meter(1.1):	Permanent display of actual lamp-current
Lamp current presetting (1.2):	Setting of lamp current in the range of 0 - 1 A.
Timer (1.3):	Shows time of device operation. Not resettable.
Power switch (1.4):	Power switch of the equipment. After turning on the device it ignites the lamp.

#### 4.2 Back Panel

**Fig. 2 Back Panel**

- Lamp (2.1):                      Lamp-connector.
- Power connector (2.2):        AC inlet with integrated fuse-holder and voltage selector

## 5 Installation

The lamp has to be connected to the supply by qualified personnel only.

**Attention!**

**Since there is high voltage up to 600 volts (even some minutes after power off) be very carefull while touching the connector.**

When setting up the instrument be sure that the ventilator and the air inlets are not covered.

Ensure that the voltage value indicated at the rear panel coincides with the voltage from your power outlet. If not, follow this procedure.

- a) Open the fuse holder and remove voltage selector.
  - b) Re-insert the voltage selector according to your voltage.
  - c) Choose the right fuse (see Technical Specifications)
- Connect the device to the power outlet with the delivered power cord.

The lamp is to be connected with the enclosed lamp cable.

## 6 Operation

After turning on the device the lamp will be ignited by an ignition voltage of approx. 600V. Then the lamp current will be regulated to the presetted value (between 0 - 1A). While operating the actual lamp-current is displayed on the panelmeter.

## 7 Pin assignment

The lamp connector has this pin assignment:

Pin	Assignment
1	not used
2	lamp, anode
3	lamp, cathode
	protective earth

## 8 Precautions

To prevent fire or shock, do not expose the unit to rain or moisture. To avoid electrical shock, do not open the cabinet. The ventilator and the air-inlets should not be covered. Refer servicing to qualified personnel only.

**NOTICE FOR COSTUMERS IN THE U.K.  
 IMPORTANT**

The wires in this mains lead are coloured in accordance with the following code:

Blue: Neutral

Black or Brown: Live

As the colours of the wires in the mains lead of this apparatus may not correspond with the coloured markings identifying the terminals in your plug proceed as follows: The wire which is coloured blue must be connected to the terminal which is marked with the letter N or coloured black. The wire which is coloured black or brown must be connected to the terminal which is marked with the letter L or coloured red.

Do not connect either wire to the earth terminal in the plug which is marked by the letter E or by the safety earth symbol or coloured green or green-and-yellow.